

LIFT BATH TROLLEY

Bolero



A reliable lift bath trolley

Bolero™ makes assisted bathing routines for dependent residents a one-to-one health care process.

For the resident, being bathed by a single caregiver is less stressful. *Bolero* provides a secure environment throughout the bathing cycle, allowing the resident to relax and enjoy the full benefits of bathing. After retrieval from the bed, the resident remains safe on the lift bath trolley during transport, transfers and the bathing session.

A single caregiver can safely manage all stages of the bathing cycle. The efficient functions and features of *Bolero* ensure a safe, ergonomic working environment and give the caregiver more time to assist the resident. With a pediatric version also available, *Bolero* can enhance the safety and comfort of bathing routines for dependent residents of all ages.

Key features



Easy transfers

The resident is easily rolled or laterally transferred onto the stretcher, which can be adjusted to the correct working height.



Bathing

Bolero integrates with Arjo recumbent bathing systems.



Folding side rail

The stretcher can be fitted with a foldable side rail to enhance resident safety and comfort.



Quick-charge battery

The quick-charge battery system ensures *Bolero* is always available for use.



Electronic scale

Bolero can also be used to weigh residents using the optional foot-operated electronic scale.



Handles

Handles on the pillar make the lift bath trolley easy to maneuver and position.

PRODUCT INFORMATION

Lifting capacity	300 lbs (136 kg)
Weight	154 lbs (70 kg)
Width	32½" (825 mm)
Length	72¾ - 76¾" (1,850-1,950 mm)
Length (pediatric)	53" - 57" (1,350-1,450 mm)
Maximum total height	67¾" (1,720 mm)
Lifting stroke range	22¼" (565 mm)
Powered functions	Raising and lowering
Bath compatibility	Rhapsody®/Primo® (System 23), Harmonie, and <i>Rhapsody</i> (System 25)
Best suited for	Mobility levels Doris and Emma



Arjo_HYG_Bolero_slick.01.0.NA.EN

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure injuries and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind. www.arjo.com

Arjo Inc. · 2349 West Lake Street · Addison, Illinois 60101 · 800-323-1245 · www.arjo.us

Arjo Canada Inc. · 90 Matheson Blvd West, Suite 350 · Mississauga, ON, L5R 3R3 · 800-665-4831 · Info.Canada@arjo.com · www.arjo.ca

arjo