

ARJOHUNTLEIGH

GETINGE GROUP

CAREVO SHOWER TROLLEY



CAREVO™ shower trolley from ArjoHuntleigh is an innovative hygiene solution for immobile patients.

The new features, innovative design, and intuitive user interface provide an easy to use and efficient hygiene solution for both caregivers and immobile patients.

Designed with Emma (E) in mind.




...with people in mind

www.ArjoHuntleigh.com


KEY FEATURES




- Accessible at all times
- Simple, intuitive and efficient



- Enhanced sense of security
- Promotes participation and mobility
- Provides support




- Improved comfort
- Pressure relieving
- Easier handling and access



- Work closer to your patient
- Reduce static load
- Improve work-flow

Product Information	
Max total load (SWL):	200 kg(440 lbs)
Max patient weight:	182 kg(401 lbs)
Weight (Carevo with mattress):	82 kg(181 lbs)
Width (stretcher):	586 mm(23")
Total width (incl. side supports):	759 mm(30")- 888 mm(35")
Total length:	2005 mm(79")
Length of chassis (incl. wheel frame):	1203 mm(47½")



5 Elevated back support
(Optional)

- Optimal patient comfort
- Easier and efficient patient transfers




6 Width adjustment option

- Greater versatility for more patient requirements
- Easier transfers and more handling options



7 Integrated shower handle holder

- Intuitive
- Easy to reach



8 Optimized head pillow and softer mattress

- Superior patient comfort
- Durable, easy to clean and disinfect



9 Simple daily maintenance

- Easy mattress attachment
- Easy mattress removal
- Easy to disinfect



Product Information

Width of chassis:	700 mm(27½")
Height of chassis (over wheels):	204 mm(8")
Height, top position floor to stretcher:	955 mm(37½")
Height, low position floor to stretcher:	605 mm(23¾")
Stroke range:	350 mm(13¾")
Drain hose holders included:	2 pieces
Max back support angle (Optional):	15°
Equipped with head pillow and foldable side supports with safety lock	

EXPLORE THE FULL RANGE OF ARJOHUNTLEIGH SHOWER SOLUTIONS AND SEE CAREVO IN ACTION

Visit www.ArjoHuntleigh.com/carevo



the preventable .com

Clinical research has shown that shower trolley solutions, like CAREVO, reduce the risk of serious injury for both caregivers and patients during hygiene tasks. For more information, visit the **preventable.com**

ARJOHUNTLEIGH

GETINGE GROUP

www.ArjoHuntleigh.co.uk

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of **ArjoHuntleigh**, **GETINGE** and **MAQUET**. **ArjoHuntleigh** focuses on patient mobility and wound management solutions. **GETINGE** provides solutions for infection control within healthcare and contamination prevention within life sciences. **MAQUET** specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.

ArjoHuntleigh UK – Head Office
ArjoHuntleigh House,
Houghton Hall Business Park,
Houghton Regis, Bedfordshire,
LU5 5XF, United Kingdom
Sales Phone: +44 (0) 1582 745700
Sales Fax: +44 (0) 1582 745745
Rental 24hr Helpline
Lo-call Phone: +44 (0) 8457 342000
Rental Fax: +44 (0) 1582 745843
Email: sales.admin@ArjoHuntleigh.com

**ArjoHuntleigh UK
Service and Technical Centre**
St Catherine Street,
Gloucester, GL1 2SL,
United Kingdom
Phone: +44 (0) 8456 114114
Service Fax: +44 (0) 1452 525207
Email: ukservice@ArjoHuntleigh.com

**ArjoHuntleigh UK
Belfast Office**
19 Heron Road,
Sydenham Business Park,
Belfast, BT3 9LE,
Northern Ireland
Phone: +44 (0) 28 9050 2000
Fax: +44 (0) 28 9050 2001
Email: belfast@ArjoHuntleigh.com
www.ArjoHuntleigh.com

ArjoHuntleigh (Ireland) Ltd.
EA House,
Damastown Industrial Park,
Mulhuddart, Dublin 15, Ireland
Phone: +353 (0)1 809 8960
Fax: +353 (0) 1 809 8971
Email: dublin@ArjoHuntleigh.com
www.ArjoHuntleigh.com