



Functionality,
meets style

ACUTE & LONG TERM CARE RANGE

Prioma™ Ward Furniture

Designed to compliment your healthcare setting



Convenient and practical

The Prioma range of bedside cabinets and overbed tables provide excellent extensions and additions to your hospital beds. The tables can be used both when the patients are in bed or seated in a chair, and provide invaluable support for everyday activities, such as reading, eating and drinking.

Blending comfort and style, with functionality and quality

The total look and feel of your room is more than just a colour palette, as is the comfort and functionality of the space being more than just the bed. At Arjo we design our products with the total space in mind, which is why we offer stylish and modern room furniture to compliment your entire setting and most importantly provide the feeling of comfort, security and style.



Overbed Tables

- Mild steel construction powder coated off-white
- Locks in compression - no release required for raising
- Modern ABS plastic base cover
- 55mm diameter castors
- Safe Working Load (SWL) 25Kg, evenly distributed
- Hettich bracket (OBT2 only)



Prioma OBT1
Fixed top table










Prioma OBT2
Split top table



Prioma OBT3
Fixed top table



Prioma OBT4
Fixed top table

Top colour options		
Prioma OBT1 Prioma OBT2	 	Prioma OBT3
	Sand White Sky Blue	  
		Prioma OBT4
		 
		Silver Grey Sky Blue

Colour and grain finish may vary slightly from swatches in this brochure.

Model	Top - design & material	Overall size	Height adjustment
Prioma OBT1	Fixed, thermo formed PVC	900 mm x 425 mm	745 mm to 1055 mm (+/- 5mm)
Prioma OBT2	Split with tilt, thermo formed PVC	700 mm x 425 mm and 195mm x 400 mm	745 mm to 1055 mm (+/- 5mm)
Prioma OBT3	Fixed, MDF laminated beech timber	900 mm x 425 mm	750 mm to 1060 mm (+/- 5mm)
Prioma OBT4	Fixed, high density polyethylene (HDPE)	700 mm x 425 mm and 195 mm x 400 mm	750 mm to 1060 mm (+/- 5mm)

Bedside Cabinets



Prioma BSC1



Prioma BSC2



Prioma BSC3



Prioma BSC4

Prioma BSC1 & BSC2

- Thermo formed PVC top, over MDF timber core
- MDF timber laminated side infill panels
- Protective aluminium corner extrusions
- Drawer and cupboard face MDF timber laminated
- Pull-out drawer (internal) injection moulded ABS
- 40 mm twin castors
- Towel hooks both sides

Prioma BSC3 & BSC4

- ABS injection moulded top and side panels
- Side panels drop down towel rail and hooks
- Pull out drawer doubles as flip over table surface
- Two brake and two free wheeling castors
- Open bottom shoe shelf (BSC4 only)

Model	Top & panel type	Overall size	Castors
Prioma BSC1	Thermo formed PVC top with MDF laminated side panels	500 mm (l) x 495 mm (w) x 870 mm (h)	Free castors, no lock
Prioma BSC2	Thermo formed PVC top with MDF laminated side panels	500 mm (l) x 495 mm (w) x 870 mm (h)	Free castors, no lock
Prioma BSC3	ABS injection moulded top and side panels	460 mm (l) x 450 mm (w) x 725 mm (h)	Free castors, with lock
Prioma BSC4	ABS injection moulded top and side panels	460 mm (l) x 450 mm (w) x 860 mm (h)	Free castors, with lock

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind. www.arjo.com

Arjo AB, Hans Michelsengatan 10, SE-211 20 Malmö, Sweden. Phone: +46 (0) 10 335 45 00

www.arjo.com

