edel immersys

Birth Pool in a Box *Professional*

Birth Pool in a Box Professional

is designed to provide a comfortable and safe environment for labouring and birthing in water for mother, baby and midwife. The single use custom fit liner acts as an infection control barrier and fits securely on the pool, which can be used multiple times and has a warranty for manufacturing faults of the lesser of 40 uses or 4 years.



Features

| Height adjustable for woman and midwife | Water depth up to 24" Integrated cupholder to help with hydration Inflation time 10 minutes | |
|--|---|--|
| 6 handles for control during labour and birth | | |
| Internal seat to sit on during labour and to sit with newborn's head safely out of the water | | |
| | Usage tracking matrix for compliance | |

| Feature | The Mini | The Regular |
|----------------------------------|---|--|
| External dimensions | 65" x 57" - 165 x 145cm | 76" 65" - 193 x 165cm |
| Internal dimensions | 45" x 37" - 114 x 94cm | 56" x 45" - 142 x 114cm |
| Internal Height | 26" - 66cm | 27" - 69cm |
| Depth of Water | Min 18" - 45cm Max 22" - 55cm | Min 18" - 45cm Max 23" - 58cm |
| Rim Width | 10" - 25cm | 10" - 25cm |
| Inflated Floor | 2" - 5cm | 2" - 5cm |
| Built-in Seat | Yes | Yes |
| Integrated Cup Holder | Yes | Yes |
| Disposable Liner | Not included - liners available separately | Not included - liners available separately |
| MyAnchor Compatible | Yes | Yes |
| Pool Material | 0.38mm Eco PVC Cadmium, Lead and Latex free | |
| Liner Material | 0.20mm Eco PVC Cadmium, | Lead and Latex free |
| Recommended Max Height of Mother | 6' or 182cm | no limit |
| Room for Partner | No | Yes |
| Volume of Water (80% Full) | 480 litres | 650 litres |
| Weight when Filled | 487kg (1,017 lbs) plus mother | 659kg (1,450 lbs) plus mother |

Materials Information

Material: PVC compliant with EU limits of phthalates Thickness: Pool 0.55mm Liner 0.20mm

Allergen Information: Latex free. Skin reaction to PVC possible but rare

Storage Temperature: min -10°C max 50°C Handling Temperature: min 20°C max 50°C

Material tested to conform to: (a) EN71 and 2005/84/EEC limiting certain phthalates, and (b) RoHS EU regulations governing cadmium, lead and other toxic substances FDA Product Code OVX, class not yet assigned.

Warranty

Covers pool for 4 years/40 uses against manufacturing fault

Related Items

Single use liners: Pro Mini Liner, Pro Regular Liner Covers: ClearFit Cover Mini, Regular Air Inflate-Deflate Pump Submersible Water Pump

